

F4/80 M-pluriBead® anti-ms

Order No.:29-F4800-20



Description

F4/80 pluriBead® suspension for the positive cell isolation of F4/80+ target cells from any murine sample material. (2ml)

Informations complémentaires

Taille	2 ml
Espèces	Souris (Murine)
Types de cellules	Monocytes, Macrophages
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.