

CD193 M-pluriBead® anti-hu

Order No.:19-19300-20



Description

CD193 pluriBead® suspension for the positive cell isolation of CD193+ target cells from any human sample material. (2ml)

Informations complémentaires

Taille	2 ml
Espèces	Humain
Types de cellules	Granulocytes basophiles, Granulocytes éosinophiles
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.