

## CD11b S-pluriBead® anti-ms

**Order No.:**29-011b0-10



### Description

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any murine sample material. (2ml)

### Informations complémentaires

Taille	2 ml
Espèces	Souris (Murine)
Types de cellules	Monocytes, Granulocytes, Macrophages, Cellules NK (sous-ensemble)
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.