

CD19 S-pluriBead® anti-ms

Order No.:29-01900-10



Description

CD19 pluriBead® suspension for the positive cell isolation of CD19+ target cells from any murine sample material. (2ml)

Informations complémentaires

Taille	2 ml
Espèces	Souris (Murine)
Types de cellules	Cellules B
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.