

## CD47 S-pluriBead® anti-sheep

**Order No.:**24-04700-10



### Description

CD47 S-pluriBead® suspension for the positive cell isolation of CD47+ target cells from any sheep sample material. (2ml)

### Informations complémentaires

Taille	2 ml
Espèces	Moutons
Types de cellules	Plaquettes, Leucocytes, Erythrocytes, RBC (globules rouges)
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.