

## CD47 S-pluriBead® anti-hu

**Order No.:**19-04700-10



### Description

CD47 pluriBead® suspension for the positive cell isolation of CD47+ target cells from any human sample material. (2ml)

### Informations complémentaires

|                            |   |
|----------------------------|---|
| Taille                     | 2 ml  |
| Espèces                    | Humain  |
| Types de cellules          | Plaquettes, Leucocytes, Erythrocytes, RBC (globules rouges)                   |
| Échantillon de matériel    | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Pureté                     | >90%  |
| Les conditions de stockage | 4-8 °C  |
| Déclaration réglementaire  | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.