

## CD117 S-pluriBead® anti-hu

**Order No.:**19-11700-10



### Description

CD117 pluriBead® suspension for the positive cell isolation of CD117+ target cells from any human sample material. (2ml)

### Informations complémentaires

Taille	2 ml
Espèces	Humain
Types de cellules	Cellules progénitrices hématopoïétiques
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.