

## CD5 S-pluriBead® anti-hu

**Order No.:**19-00500-10



### Description

CD5 pluriBead® suspension for the positive cell isolation of CD5+ target cells from any human sample material. (2ml)

### Informations complémentaires

Taille	2 ml
Espèces	Humain
Types de cellules	Les cellules T, Cellules B (sous-ensemble), Thymocytes
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.