

CD7 S-pluriBead® anti-hu

Order No.:19-00700-10



Description

CD7 pluriBead® suspension for the positive cell isolation of CD7+ target cells from any human sample material. (2ml)

Informations complémentaires

| | |
|----------------------------|---|
| Taille | 2 ml |
| Espèces | Humain |
| Types de cellules | Cellules NK, Les cellules T, Thymocytes |
| Échantillon de matériel | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Pureté | >90% |
| Les conditions de stockage | 4-8 °C |
| Déclaration réglementaire | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.