

CD8 S-pluriBead® anti-hu

Order No.:19-00800-10



Description

CD8 pluriBead® suspension for the positive cell isolation of CD8+ target cells from any human sample material. (2ml)

Informations complémentaires

Taille	2 ml
Espèces	Humain
Types de cellules	Cellules NK (sous-ensemble), Cellules T cytotoxiques
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.