

## 10% BSA Stock Solution in 1x PBS






Order No.	Size	Sterility
60-00020-10	50 ml	Stérile
60-00020-11	250 ml	Stérile
60-00020-12	500 ml	Stérile

### Description

10% BSA stock solution in 1x PBS can used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

### Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

### Informations complémentaires

Délai de livraison	1-2
Domaine d'intérêt	Biologie cellulaire
Stérilité	Stérile
Gamme de ph	pH 7.4
Conditions d'expédition	Température ambiante
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

