

## pluriSpin® Human CD45<sup>+</sup> Cell Depletion (10 ml)

**Order No.:**19-09002-32



### Description

pluriSpin® Human CD45<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched CD45<sup>+</sup>-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45<sup>+</sup> cells, while the isolated population free of CD45<sup>+</sup> are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45<sup>+</sup> cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Informations complémentaires

Taille	10 ml
Espèces	Humain
Types de cellules	CD45- Cellule
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.