

pluriSpin® Human CD4⁺ Cell Depletion



Order No.	Size
19-03002-31	2 ml
19-03002-32	10 ml

Description

pluriSpin® Human CD4⁺ Cell Depletion is designed to isolate pure, viable and untouched helper t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Helper T cells, while the isolated human helper t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified helper t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Informations complémentaires

Espèces	Humain
Types de cellules	Cellules T auxiliaires
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.