

pluriSpin® Human Monocyte Depletion


Order No.
Size

19-01002-31

2 ml

19-01002-32

10 ml

Description

pluriSpin® Human Monocyte Depletion is designed to isolate pure, viable and untouched monocytes-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels monocytes, while the isolated human monocytes-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Informations complémentaires

Espèces	Humain
Types de cellules	Monocytes
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.