

## pluriSpin® Human PBMC Depletion (10 ml)

**Order No.:**19-08002-32



### Description

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

### Informations complémentaires

Taille	10 ml
Espèces	Humain
Types de cellules	Cellules non-PBMC
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.