

## Bundle pluriStrainer® Mini 40 µm and Flow Cytometry Tubes (25 pcs.)

**Order No.:**43-10040-46



### Description

Bundle includes sterile pluriStrainers Mini with 40 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 40 µm polyethylenterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 40 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

### Candidature

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marcol-sm-12, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry(FACS™)

### Informations complémentaires

Délai de livraison	1-2
Taille des mailles	40 µm
Emballage	25 pcs. (25 pcs. / Bag)
Taille	25 pcs.
Couleur	bleu clair
Contrôle des flux	Non
Mesh Fixation	Injected in housing for strong hold
Tissu de la matière	PET (polyéthylène téréphthalate)

Matériau du boîtier	PP (Polypropylène)
S'inscrit dans	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
Traitement	Gripping surface enables easy and aseptic handling.
Stérilité	Stérile
Comparable à	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
Stabilité	RLT® Buffer, TRIzol®, isopropanol, organic solvents
Centrifugable	Oui
Dissociation des tissus	Oui
Désinfectable	Oui
Conditions d'expédition	Température ambiante
Les conditions de stockage	Température ambiante
Informations juridiques	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.