

Centrifuge Tube 50 ml (sterile, 300 pcs.)



Order No.:05-00001-01

Description

50 ml tubes are suitable for a wide range of laboratory applications such as sample preparation or sample storage. The tubes are ideal for the use with our **pluriStrainers** and guarantee a perfect fit with the Connector Ring. Moreover the tubes can be used as container for **density gradient centrifugation** or the enrichment of primary cells with **pluriBead**.

Screw tube, 50 ml, 114 x 28 mm, pointed bottom, PP, with printed writing field and scaling, with mounted red cap, sterile and pyrogen-free / endotoxin-free, 25 pcs./bag

Candidature

- Mixing
- Incubation
- Centrifugation
- Storage
- Collection vessel for filtration and purification of liquids

Informations complémentaires

Taille	300 pcs.
Pré-rempli	Unfilled
Délai de livraison	1-2
Variation	20 pcs, sterile
Emballage	25 pcs. / Bag
Les conditions de stockage	Température ambiante
Conditions d'expédition	Température ambiante

Déclaration réglementaire

For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.