

Funnel (50 pcs. , non-sterile)



Order No.:42-50000-50

Description

Funnel - Sterile device to support sample load onto pluriStrainer.

Informations complémentaires

Taille	50 pcs.
Délai de livraison	1-2
Variation	50 pcs, non-sterile
Conditions d'expédition	Température ambiante
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.