

ÜberStrainer Set 2

Order No.:43-70000-98



Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 1 µm, 5 µm, 10 µm, 30 µm, 60 µm.

Informations complémentaires

Taille des mailles	Multiple
Taille	10 pcs
Tissu de la matière	PET (polyéthylène téréphtalate)
Matériau du boîtier	PP (Polypropylène)
Mesh Type	woven
Stérilité	Stérile
Conditions d'expédition	Température ambiante
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.