

Bundle pluriStrainer® Mini 100 µm and Flow Cytometry Tubes (100 pcs.)

Order No.:43-10100-66



Description

Bundle includes sterile pluriStrainers Mini with 100 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 100 µm polyethylenterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 100 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

Candidature

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marrow, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry(FACS™)

Informations complémentaires

| | |
|---------------------|-------------------------------------|
| Délai de livraison | 1-2 |
| Taille des mailles | 100 µm |
| Emballage | 100 pcs. (25 pcs. / Bag) |
| Taille | 100 pcs. |
| Couleur | jaune |
| Contrôle des flux | Non |
| Mesh Fixation | Injected in housing for strong hold |
| Tissu de la matière | PET (polyéthylène téréphthalate) |

| | |
|----------------------------|---|
| Matériau du boîtier | PP (Polypropylène) |
| S'inscrit dans | 1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate |
| Traitement | Gripping surface enables easy and aseptic handling. |
| Stérilité | Stérile |
| Comparable à | Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters |
| Stabilité | RLT® Buffer, TRIzol®, isopropanol, organic solvents |
| Centrifugable | Oui |
| Dissociation des tissus | Oui |
| Désinfectable | Oui |
| Conditions d'expédition | Température ambiante |
| Les conditions de stockage | Température ambiante |
| Informations juridiques | Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.