

## Syringe Strainer 85 µm (12 pcs., non-sterile)

**Order No.:**43-71085-50



### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Informations complémentaires

|                            |  |
|----------------------------|--|
| Taille des mailles         | 85 µm  |
| Variation                  | 12 pcs, non-sterile  |
| Taille                     | 12 pcs.  |
| Contrôle des flux          | Oui  |
| Matériau du boîtier        | LD-PE (polyéthylène basse densité)                           |
| Mesh Type                  | woven  |
| Stérilité                  | Stérile  |
| Conditions d'expédition    | Température ambiante   |
| Les conditions de stockage | Température ambiante   |
| Déclaration réglementaire  | For research use only. Not for use in diagnostic procedures. |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.