

Syringe Strainer 70 µm (12 pcs., sterile)

Order No.:43-71070-51



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Informations complémentaires

Taille	12 pcs.
Contrôle des flux	Oui
Taille des mailles	70 µm
Matériau du boîtier	LD-PE (polyéthylène basse densité)
Stérilité	Stérile
Variation	12 pcs. sterile
Conditions d'expédition	Température ambiante
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

