

## Syringe Strainer 20 µm





Order No.	Sterility	Variant
43-71020-50	Stérile	12 pcs, non-sterile
43-71020-51	Stérile	12 pcs sterile

### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Informations complémentaires

Taille des mailles	20 µm
Taille	12 pcs.
Contrôle des flux	Oui
Matériau du boîtier	LD-PE (polyéthylène basse densité)
Mesh Type	woven
Stérilité	Stérile
Conditions d'expédition	Température ambiante
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

