

Syringe Strainer 300 µm



Order No.	Sterility	Variant
43-71300-50	Stérile	12 pcs, non-sterile
43-71300-51	Stérile	12 pcs sterile

Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Informations complémentaires

Taille des mailles	300 µm
Taille	12 pcs.
Contrôle des flux	Oui
Matériau du boîtier	LD-PE (polyéthylène basse densité)
Mesh Type	woven
Stérilité	Stérile
Conditions d'expédition	Température ambiante
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

