

pluriSpin® Human T Cell Enrichment (2 ml)

Order No.:19-02001-31



Description

pluriSpin® Human T Cell Enrichment is designed to isolate pure, viable and untouched T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, B cells, NK cells and granulocytes, while the negatively isolated human T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified T cells are present as a highly enriched population at the interface between the plasma and the density medium.

Informations complémentaires

Taille	2 ml
Espèces	Humain
Types de cellules	Les cellules T
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

