

Re-Strainer (In-Line-Strainer) 1 µm - Sterile

Order No.:43-75001-40



Description

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

Candidature

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

Informations complémentaires

Délai de livraison	1-2
Taille des mailles	1 µm
Variation	20 pcs, sterile
Emballage	4 devices per bag

Taille	20 pcs.
Couleur	rouge
Contrôle des flux	oui
Mesh Fixation	Injected in housing for strong hold
Tissu de la matière	PET (polyéthylène téréphtalate)
Matériau du boîtier	LD-PE (polyéthylène basse densité)
Stérilité	Stérile
Stabilité	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
Conditions d'expédition	Room Temperature
Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only
Informations juridiques	The product is for research and development only, not for diagnostic or therapeutic use

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.