

pluriSpin® Human Monocyte Depletion (2 ml)

Order No.:19-01002-31



Description

pluriSpin® Human Monocyte Depletion is designed to isolate pure, viable and untouched monocytes-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels monocytes, while the isolated human monocytes-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Informations complémentaires

| | |
|----------------------------|---|
| Taille | 2 ml |
| Espèces | Humain |
| Types de cellules | Monocytes |
| Échantillon de matériel | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Pureté | >90% |
| Les conditions de stockage | 4-8 °C |
| Déclaration réglementaire | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

