

## 10x PBS Solution (500 ml)

**Order No.:**60-00010-11



### Description

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- RNase/DNase free
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (10mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

大小	500 ml
不育症	Sterile
相位范围	pH 6.8
运输条件	Room Temperature
储存条件	Room Temperature

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.