

EDTA Solution (50 ml)

Order No.: 60-00030-11



Description

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as an chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

Additional Information

大小	50 ml
感兴趣的领域	Cell Biology
不育症	Sterile
相位范围	pH 7.4
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.