

pluriPlix



Order No.:50-01010-80

Description

Universal rotation mixer without own drive system for usage on top of a magnetic stirrer (not included in delivery). Particularly designed for the gentle mixing of cells. Allows for 2 mixing options: over-head and tilting-rotation (plus: adjustable rotation angle).

Additional Information

大小	1 device
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.