

Leuko Spin Medium (250 ml)

Order No.:60-00091-11



Description

Leuko Spin is a sterile, ready-to-use density medium for the isolation of all Leukocytes from fresh human peripheral blood or buffy coat with a high yield. Mononuclear and polymorphonuclear (multi-lobed nuclei) cells can be enriched by a single-step density gradient centrifugation with Leuko Spin Medium. Through combination with PBMC Spin in a double gradient centrifugation it is possible to isolate and separate granulocytes and PBMCs in one step with high yield and purity. The enriched cells can be used for a wide range of downstream applications such as MACS (magnetic activated cell sorting, a trademark of Miltenyi Biotec GmbH) or stem cell research.

Additional Information

| | |
|-------|--|
| 大小 | 250 ml |
| 交付时间 | 1-2 |
| 储存条件 | Room temperature, protect from light |
| 规范性声明 | For research use only. Not for use in diagnostic procedures. |
| 相位范围 | 6.8 - 7.4 |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.