

CD4 M-pluriBead® anti-ms

Order No.:29-00400-20



Description

CD4 pluriBead® suspension for the positive cell isolation of CD4+ target cells from any murine sample material. (2ml)

Additional Information

大小	2 ml
品种	Mouse (Murine)
细胞类型	Monocytes (subset), Helper T cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.