

CD8 M-pluriBead® anti-ms

Order No.:29-00800-20



Description

CD8 pluriBead® suspension for the positive cell isolation of CD8+ target cells from any murine sample material. (2ml)

Additional Information

大小	2 ml
品种	Mouse (Murine)
细胞类型	NK cells (subset), Cytotoxic T cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.