

## CD117 M-pluriBead® anti-hu

**Order No.:**19-11700-20



### Description

CD117 pluriBead® suspension for the positive cell isolation of CD117+ target cells from any human sample material. (2ml)

### Additional Information

大小	2 ml
品种	Human
细胞类型	Hematopoietic progenitor cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.