

CD7 M-pluriBead® anti-hu

Order No.:19-00700-20



Description

CD7 pluriBead® suspension for the positive cell isolation of CD7+ target cells from any human sample material. (2ml)

Additional Information

大小	2 ml
品种	Human
细胞类型	NK cells, T cells, Thymocytes
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.