

PLT Spin Medium (250 ml)

Order No.:60-00094-11



Description

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

Additional Information

大小	250 ml
交付时间	1-2
储存条件	Room temperature, protect from light
规范性声明	For research use only. Not for use in diagnostic procedures.
相位范围	6.8 - 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.