

## CD31 S-pluriBead® anti-ms

**Order No.:**29-03100-10



### Description

CD31 pluriBead® suspension for the positive cell isolation of CD31+ target cells from any murine sample material. (2ml)

### Additional Information

大小	2 ml
品种	Mouse (Murine)
细胞类型	Platelets, Endothelial cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.