

CD55 S-pluriBead® anti-hu

Order No.:19-05500-10



Description

CD55 pluriBead® suspension for the positive cell isolation of CD55+ target cells from any human sample material. (2ml)

Additional Information

大小	2 ml
品种	Human
细胞类型	Hematopoietic cells, Non-hematopoietic cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.