

CD271 S-pluriBead® anti-hu

Order No.:19-27100-10



Description

CD271 pluriBead® suspension for the positive cell isolation of CD271+ target cells from any human sample material. (2ml)

Additional Information

大小	2 ml
品种	Human
细胞类型	Mesemchymal cells, Schwann cells, Stem/Stromal cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.