

## 10% BSA Stock Solution in 1x PBS






Order No.	Size	Sterility
60-00020-10	50 ml	Sterile
60-00020-11	250 ml	Sterile
60-00020-12	500 ml	Sterile

### Description

10% BSA stock solution in 1x PBS can used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

### Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

### Additional Information

交付时间	1-2
感兴趣的领域	Cell Biology
不育症	Sterile
相位范围	pH 7.4
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

