

## EDTA Solution (0.2M in PBS)






Order No.	Size	Sterility
60-00030-10	10 ml	Sterile
60-00030-11	50 ml	Sterile
60-00030-12	100 ml	Sterile

### Description

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as a chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

### Additional Information

感兴趣的领域	Cell Biology
不育症	Sterile
相位范围	pH 7.4
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.