

EDTA Solution (0.2M in PBS)


☐
☐
☐

Order No.	Size	Sterility
60-00030-10	10 ml	Sterile
60-00030-11	50 ml	Sterile
60-00030-12	100 ml	Sterile

Description

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as an chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

Additional Information

感兴趣的领域	Cell Biology
不育症	Sterile
相位范围	pH 7.4
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.