

## pluriSpin® Human CD4<sup>+</sup> Cell Depletion





Order No.	Size
19-03002-31	2 ml
19-03002-32	10 ml

### Description

pluriSpin® Human CD4<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched helper t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Helper T cells, while the isolated human helper t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified helper t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

品种	Human
细胞类型	Helper T cells
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.