

pluriSpin® Human CD45⁺ Cell Depletion


Order No.
Size

19-09002-31

2 ml

19-09002-32

10 ml

Description

pluriSpin® Human CD45⁺ Cell Depletion is designed to isolate pure, viable and untouched CD45⁺-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45⁺ cells, while the isolated population free of CD45⁺ are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45⁺ cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

品种	Human
细胞类型	CD45- Cell
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.