

pluriSpin® Human PBMC Depletion



Order No.	Size
19-08002-31	2 ml
19-08002-32	10 ml

Description

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

Additional Information

品种	Human
细胞类型	Non-PBMC cells
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.