


PLT Spin Medium



Order No.

Size

Order No.	Size
60-00094-10	100 ml
60-00094-11	250 ml
60-00094-12	500 ml

Description

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

Application

Centrifugation scheme for a platelet preparation with PLT Spin Medium




Fig. 1 Preparation of Platelets with PLT Spin density separation medium.




Fig. 2 Enriched platelets are found on the upper layer of the PLT Spin Medium.

Additional Information

交付时间	1-2
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
储存条件	Room temperature, protect from light
规范性声明	For research use only. Not for use in diagnostic procedures.
相位范围	6.8 - 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.