

## ÜberStrainer 5 µm



**Order No.:**43-70005-03

### Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

### Additional Information

交付时间	30-32
Mesh Size	5 µm
变化	25 pcs, sterile, PET-Mesh
大小	10 pcs
Fabric material	PET (Polyethylenterephthalat)
外壳材料	PP (Polypropylene)
Mesh Type	woven
不育症	Sterile
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.,

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.