

## Bundle pluriStrainer® Mini 10 µm and Flow Cytometry Tubes (500 pcs.)

**Order No.:**43-10010-76



### Description

Bundle includes sterile pluriStrainers Mini with 10 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 10 µm polyethyleneterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 10 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

### Application

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marcol-sm-12, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry(FACS™)

### Additional Information

交付时间	1-2
Mesh Size	10 µm
包装	500 pcs. (25 pcs. / Bag)
大小	500 pcs.
颜色	turquoise
流量控制	No
Mesh Fixation	Injected in housing for strong hold
Fabric material	PET (Polyethyleneterephthalat)

外壳材料	PP (Polypropylene)
符合	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
处理	Gripping surface enables easy and aseptic handling.
不育症	Sterile
相当于	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
稳定性	RLT® Buffer, TRIzol®, isopropanol, organic solvents
离心	Yes
组织解离	Yes
可消毒	Yes
运输条件	Room Temperature
储存条件	Room Temperature
法律信息	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.