

## Bundle pluriStrainer® Mini 40 µm and Flow Cytometry Tubes (100 pcs.)

**Order No.:**43-10040-66



### Description

Bundle includes sterile pluriStrainers Mini with 40 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 40 µm polyethyleneterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 40 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

### Application

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marcol-sm-12, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry(FACS™)

### Additional Information

交付时间	1-2
Mesh Size	40 µm
包装	100 pcs. (25 pcs. / Bag)
大小	100 pcs.
颜色	light blue
流量控制	No
Mesh Fixation	Injected in housing for strong hold
Fabric material	PET (Polyethyleneterephthalat)

**pluriSelect Life Science**

Deutscher Platz 5a, 04103 Leipzig, DE

Phone: +49 341 333 858 0 | Fax: +49 341 333 858 99

support@pluriselect.com

<https://www.pluriselect.com/>

外壳材料	PP (Polypropylene)
符合	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
处理	Gripping surface enables easy and aseptic handling.
不育症	Sterile
相当于	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
稳定性	RLT® Buffer, TRIzol®, isopropanol, organic solvents
离心	Yes
组织解离	Yes
可消毒	Yes
运输条件	Room Temperature
储存条件	Room Temperature
法律信息	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.