

## Syringe Strainer 20 µm (12 pcs., sterile)

**Order No.:**43-71020-51



### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Mesh Size	20 µm
变化	12 pcs sterile
大小	12 pcs.
流量控制	Yes
外壳材料	LD-PE (Low Density Polyethylen)
Mesh Type	woven
不育症	Sterile
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.